

**UNCLASSIFIED**  
**AD 419021**

**DEFENSE DOCUMENTATION CENTER**

**FOR**

**SCIENTIFIC AND TECHNICAL INFORMATION**

**CAMERON STATION, ALEXANDRIA, VIRGINIA**



**UNCLASSIFIED**

**NOTICE:** When government or other drawings, specifications or other data are used for any purpose other than in connection with a definitely related government procurement operation, the U. S. Government thereby incurs no responsibility, nor any obligation whatsoever; and the fact that the Government may have formulated, furnished, or in any way supplied the said drawings, specifications, or other data is not to be regarded by implication or otherwise as in any manner licensing the holder or any other person or corporation, or conveying any rights or permission to manufacture, use or sell any patented invention that may in any way be related thereto.

64-5

NOLC REPORT 589  
15 AUGUST 1963

**AN INSTRUMENT FOR MEASURING  
THE SPECTRAL EMITTANCE OF A SOLID  
IN THE INFRARED FROM 16 TO 45 MICRONS  
AT LOW TEMPERATURES**

**D. L. STIERWALT**

**D. D. KIRK**

**J. B. BERNSTEIN**

**RESEARCH DEPARTMENT**

DDC AVAILABILITY NOTICE

Qualified requesters may obtain copies of  
this report from DDC



**NAVAL ORDNANCE LABORATORY CORONA**  
CORONA, CALIFORNIA

419021  
419021  
419021

# **NAVAL ORDNANCE LABORATORY CORONA**

**W. R. KURTZ, CAPT., USN**  
Commanding Officer

**F. S. ATCHISON, Ph. D.**  
Technical Director

## **FOREWORD**

The work described in this report was conducted in the Infrared Division, Research Department, Naval Ordnance Laboratory Corona, with funds supplied by the National Aeronautics and Space Administration Purchase Order No. W-11, 400-B.

**CURTIS J. HUMPHREYS**  
Head, Research Department

## **ABSTRACT**

As part of a continuing research project for measuring the emittance of solids, a Beckman IR-3 spectrophotometer has been modified. This modification permits the measurement of spectral emittance of solids in the 16 to 45 $\mu$  region in the temperature range between 77 and 473°K.

## INTRODUCTION

The theory, techniques, and experimental apparatus used in acquiring the spectral emittance of opaque and transparent solids in the 2 to 15  $\mu$  region and in the temperature range between 313 and 453°K have been given in a previous report,<sup>1</sup> which described in detail the modification of a Beckman IR-3 spectrophotometer and the associated components that permitted these measurements.

Other modifications have permitted the measurement of spectral emittance in the 1.8 to 25  $\mu$  region at temperatures from 313 to 473°K.

This report describes a second modification of the instrument, which permits spectral emittance measurements in the 16 to 45  $\mu$  region and in the temperature range between 77 and 473°K.

## MODIFIED EQUIPMENT

Figure 1 shows a diagram of the optical system of the unmodified IR-3. The gas cell and liquid cell were not required for this experiment and were removed; this eliminated lenses G and L and windows S and V. Figure 2 shows the optical system after the modification described in this report. In the monochromator compartment, the Littrow prism-mirror system (M and N on Figure 1) was replaced by a Bausch and Lomb replica grating; this grating, with 20 grooves per mm, has a 64 x 64 mm ruled area, and is blazed at 45  $\mu$ . The full travel of the wavelength drive permits the grating to rotate through an angle of 18 deg; this rotation allows a spectral range from 16 to 45  $\mu$  in the first order. The increased wavelength range covered by this modification required that the wavelength drive-gear train be reduced by a factor of 4. The other Littrow prism and mirror system (M' and N' on Figure 1) was replaced by a plane mirror. Lens K was replaced with a cesium iodide window; this is possible because the optical path is shortened by the elimination of the liquid and gas cell compartments. Lens Q was removed but not replaced with a window. The curved slits were replaced with straight slits that open to a maximum of 7 mm. The

---

<sup>1</sup>NAVWEPS Report 7160, Infrared Spectral Emissivity of Optical Materials, by D. L. Stierwalt, 15 January 1961. (NOLC 537)




FIGURE 1. Unmodified Optical System




FIGURE 2. Optical System After Second Modification

selection of straight slits was made because in this energy-limited region where wide slits are necessary for obtaining relatively low resolution, little advantage can be gained from a curved slit.

The monochromator was calibrated by the use of a  $6-6.5\mu$  bandpass filter and the  $6.238\mu$  absorption line of polystyrene in the third to seventh orders. Since the calibration curve is nearly a straight line, these five points are sufficient to ensure accuracy of approximately 0.1 percent.

The source compartment was completely redesigned to permit the sample holder with sample to be inserted into the compartment at the end of a cold finger. Figure 3 is a cross-sectional view of the cold finger and sample holder. To keep thermal conduction low, the wall of the cold finger is stainless steel, 0.010 in. thick. Because the source compartment is evacuated, there is no need for a vacuum jacket around the cold finger. The sample holder with sample is inserted into the compartment in such a manner that one of the exposed sample faces is perpendicular to, and along the optical axis of, a concave mirror. This mirror focuses the energy of the sample on the entrance slit of the monochromator. The other exposed face of the sample (S on Figure 2) looks into a wedge-shaped blackbody (B on Figures 2 and 3) which is maintained at monochromator temperature. This blackbody provides a known background for transparent samples and also eliminates any radiation which might be reflected back through the sample. Other components of




FIGURE 3. Cold Finger and Sample Holder, Cross Section



the compartment are a 10-cycle chopper (C on Figure 2), a reference cavity, and two bandpass filters.

One bandpass interference filter covers the region of 14 to 25  $\mu$ ; the other, 25 to 45  $\mu$ . They can be manually adjusted from outside the source compartment at appropriate times during a measurement, to provide the proper order separation. Figure 4 shows the transmission characteristics of the filters.

The chopper has a mirror surface so that the radiation from the reference cavity enters the monochromator when the chopper interrupts the sample beam. If the reference cavity, background wedge, and sample chamber are kept at the same temperature as the monochromator, then only radiation emitted by the sample contributes to the ac signal at the detector.

For making transmission measurements, a Nernst glower can be used in place of the sample holder; the sample is then placed in front of the exit slit in the receiver compartment. The source compartment can be evacuated separately from the monochromator compartment.

Figure 5 shows the interior of the source compartment; Figure 6, the source compartment with the cold finger and sample holder in position for a measurement; and Figure 7, the cold finger with the blackbody in position in the sample holder.

In the receiver compartment, the photodetector and lead sulfide detectors were removed, and the thermocouple with a potassium bromide lens was replaced by a thermocouple with a cesium iodide lens.

The receiver and source compartments were joined directly to the monochromator compartment.

## MEASUREMENT PROCEDURE

The previous report developed the general theory and presented a comparatively intricate method for reducing the acquired data. By utilizing a special case developed from the general theory, the method for reducing the data has been simplified. In most instances, a continuous record of the emittance can be obtained rather than a point-by-point result.

From the general theory, it was shown that in the special case where the blackbody temperature is the same as the sample temperature, the





FIGURE 4. Transmission Characteristics of Bandpass Filters



FIGURE 5. Source Compartment Interior



FIGURE 6. Source Compartment With Cold Finger and Sample Holder



**FIGURE 7. Cold Finger With Blackbody Inserted in Sample Holder**

emittance of the sample is the ratio of the signal from the sample to the signal from the blackbody, or

$$e_s = \frac{S_s}{S_{bb}}$$

where

$e_s$  = emittance of the sample

$S_s$  = signal from the sample

$S_{bb}$  = signal from the blackbody

This feature of the general theory can be applied to the instrumentation in the following manner: A blackbody is inserted into the sample holder and put in the source compartment in the same position as that of a sample. A standardization of this blackbody energy, over the spectral range and at the same temperature at which the sample's emittance is to be measured, is recorded on magnetic tape and on a chart recorder. The instrumentation is so devised that, during this standardization, a program of wavelength position as a function of slit width is recorded on the tape; this will permit examination of a constant energy to the detector in the spectral region. Where the energy from the blackbody is sufficient, the spectral region can be represented on the chart as a full scale deflection and designated 100 percent; when the sample energy is being tested, the recording on the chart is a direct, continuous reading of emittance.

In the regions where the energy from the blackbody is insufficient to give a full scale deflection, energy to the detector is not constant. The ratio of the energy obtained from the sample to that obtained from the blackbody can be read from the chart, thus giving a point-by-point measurement of the sample's emittance.

Some of the first measurements made with the modified instrumentation were on samples of cadmium sulfide and silicon at 77, 203, and 473°K. Figure 8 shows the spectral emittance of a sample of CdS approximately 4 mm thick at 77 and 203°K. Figure 9 shows the emittance of a sample of Si 2 mm thick measured at 77°K.

## DISCUSSION

The purpose of this report is to show the further modifications to the Beckman IR-3 spectrophotometer that were necessary to extend the




FIGURE 8. Spectral Emittance of a Cadmium Sulfide Sample 4 mm Thick




FIGURE 9. Spectral Emittance at 77°K of a Silicon Sample 2 mm Thick



spectral region and temperature range of the instrumentation discussed in the earlier report. Over this spectral region, CdS and Si exhibit high and low emittances, respectively; data from both are included to show the new capability of the instrument.

All the solids that were measured by the instrument before this modification will be measured again in the new spectral region that has been made available. The results of these measurements will be the subjects of future reports.

NOLC Report 589, NASA Purchase Order No. W-11, 400-B.

# INITIAL DISTRIBUTION

	<u>Copies</u>		<u>Copies</u>
Chief, Bureau of Naval Weapons		Commander	
Navy Department		Naval Ordnance Test Station	
Washington 25, D. C.		China Lake, Calif.	
Attn: Code R-3	1	Attn: H. Bennett	1
R-5	1	L. Nichols	1
R-12	1	A. L. Olson	1
RM	1	B. Seraphin, Code 5019	1
RMMO	1	Technical Library	1
RMWC	1		
RRMA-3	1	Commanding Officer	
RRRE	1	Naval Ordnance Test Station	
DLI-31	4	Pasadena Annex	
		3202 E. Foothill Blvd.	
Chief, Bureau of Ships		Pasadena 8, Calif.	
Navy Department		Attn: Library	1
Washington 25, D. C.			
Attn: Code 452F	1	Commander	1
681A	1	Naval Weapons Laboratory	
687C1	1	Dahlgren, Va.	
694E, D. Wilson	1	Director	
Commander		National Security Agency	
Naval Ordnance Laboratory		Fort George G. Meade, Md.	
White Oak		Attn: REMP-2	1
Silver Spring, Md.		Superintendent	
Attn: Dr. J. Dixon, Solid State		Naval Postgraduate School	
Division	1	Monterey, Calif.	
		Attn: Library	1
Chief of Naval Research		E. C. Crittenden, Jr.	1
Navy Department			
Washington 25, D. C.		Commanding Officer & Director	
Attn: Code 407	1	Naval Radiological Defense	
420	1	Laboratory	
Commanding Officer	1	San Francisco 24, Calif.	
Office of Naval Research		Attn: L. D. Miller, Code 944A	1
Pasadena Branch			
1030 East Green Street		Commander	
Pasadena 1, Calif.		Pacific Missile Range	
		Naval Missile Center	
Director		Point Mugu, Calif.	
Naval Research Laboratory		Attn: Technical Library	1
Washington 25, D. C.			
Attn: R. Wallis	1		

	<u>Copies</u>		<u>Copies</u>
Director		Aerojet-General Corporation	
Navy Electronics Laboratory		Advanced Research Division	
San Diego, Calif.		Azusa, Calif.	
Attn: H. E. Lee	1	Attn: R. C. Carlston, Bldg. 415	1
Technical Library	1		
Officer in Charge		Avco Corporation	
Naval Civil Engineering Laboratory		Research and Development Division	
Port Hueneme, Calif.		Physics Department - T430	
Attn: Technical Library	1	Wilmington, Mass.	
		Attn: Dr. T. Wentink, Jr.	1
Officer in Charge			
Naval Weapons Services Office		Barnes Engineering Company	
U. S. Naval Station		Stamford, Conn.	
(Washington Navy Yard Annex)		Attn: R. DeWaard	1
Washington 25, D. C.			
Attn: Code 120	2	Bell and Howell Research Center	
		360 Sierra Madre Villa	
Director		Pasadena, Calif.	
National Bureau of Standards		Attn: Dr. R. K. Willardson,	
Washington 25, D. C.		Solid State	1
Attn: Dr. H. P. R. Frederikse,			
Code 404	1	Bulova Research and Development	
Dr. A. Maki	1	Laboratories	
Dr. R. D. Huntoon, Code 30	1	6210 Woodside Avenue	
		Woodside 77, Long Island, N. Y.	
Commanding General		Attn: Dr. T. K. Steele	1
Army Electronics Proving Ground			
Fort Huachuca, Ariz.		California Institute of Technology	
Attn: SIGPG-DCGES-2	1	1201 East California Street	
		Pasadena 4, Calif.	
Commander		Attn: I. E. Newland, Reports	
U. S. Army Materiel Command		Group JPL	1
Harry Diamond Laboratories			
Washington 25, D. C.		Eastman Kodak Company	
Attn: Dr. W. K. Saunders	1	Naval Ordnance Division	
		50 West Main Street	
Commander		Rochester 14, Minn.	
Cambridge Research Center		Attn: W. McKusick	1
Laurence G. Hanscom Field			
Bedford, Mass.		Ford Motor Company	
Attn: CROTLA	1	Aeronutronics Division	
J. Howard	1	Research Operations	
		Newport Beach, Calif.	
Defense Documentation Center	20	Attn: L. M. Lambert	1
Bldg. #5, Cameron Station			
Alexandria, Va.		General Electric Company	
		One River Road	
		Schenectady 5, N. Y.	
		Attn: G. C. Dodson	1

	<u>Copies</u>		<u>Copies</u>
General Telephone and Electronics Laboratory	1	Marshall Space Flight Center Huntsville, Ala.	
1015 Corporation Way		Attn: G. Heller	1
Palo Alto, Calif.		Massachusetts Institute of Technology	
Hughes Research Laboratory		Lincoln Laboratory	
Hughes Aircraft Company		Lexington 73, Mass.	
Malibu, Calif.		Attn: J. H. Chisholm	1
Attn: C. Asawa	1	National Aeronautics and Space Administration	
Illinois Institute of Technology		Goodard Space Flight Center	
Technology Center		Greenbelt, Md.	
10 West 35th Street		Attn: Mr. Catoe	1
Chicago 18, Ill.		Dr. L. Dunkelman	1
Attn: Dr. J. W. Buttrey	1	Dr. R. Hanel	1
Dr. S. Nudelman	1	Dr. W. Nordberg	1
Jet Propulsion Laboratory		National Aeronautics and Space Administration	
4800 Oak Grove Drive		Washington 25, D. C.	
Pasadena, Calif.		Attn: C. P. Mook, Code RV-1	1
Attn: R. H. Chase	1	Dr. W. A. Menzel	1
Mr. Friedman	1	Pennsylvania State College	
Mr. Powell	1	State College, Pa.	
W. Schade	1	Attn: Director, Ordnance Research Laboratory	1
The Johns Hopkins University		Royal Aircraft Establishment	
Applied Physics Laboratory		Radio Department	
8621 Georgia Avenue		Farnborough, Hants, England	
Silver Spring, Md.		Attn: T. S. Moss	1
Attn: Dr. A. I. Mahan	1	Via: BuWeps, Code DSC-3	
A. Nagy	1	Royal Radar Equipment	
Lewis Research Center		Ministry of Aviation	
Cleveland, Ohio		St. Andrews Road	
Attn: J. L. Pollack	1	Malvern	
Lockheed Aircraft Corporation		Worchestershire, England	
Plant B-1		Attn: E. H. Putley	1
Burbank, Calif.		Via: BuWeps, Code DSC-3	
Attn: Dr. L. G. Mundie,		Santa Barbara Research Center	
Dept. 7210, Bldg. 167A	1	Santa Barbara Airport	
Lockheed Aircraft Corporation		Goleta, Calif.	
Missile Systems Division		Attn: F. Young	1
Sunnyvale, Calif.			
Attn: M. Eisenberg	1		
Dr. C. Kooi	1		

	<u>Copies</u>		<u>Copies</u>
Servo Corporation of America		NOLC:	
111 New South Road		J. A. Hart, Code RT-3	1
Hicksville, Long Island, N. Y.		C. J. Humphreys, Code RR	2
Attn: I. Melman	1	R. L. Conger, Code RRM	1
Signals Research and Development		R. W. Sharp, Code RRMM	1
Establishment		R. F. Potter, Code RRI	1
Ministry of Aviation		J. B. Bernstein, Code RRIO	1
Christchurch		D. D. Kirk, Code RRIO	1
Hampshire, England		V. A. Patton, Code RRIO	1
Attn: Dr. M. R. Brown,		D. L. Stierwalt, Code RRIO	1
(ZSP/5/08)	1	G. Zaeschmar, Code RRIO	1
Via: BuWeps, Code DSC-3		R. A. Buchanan, Code RRIS	1
Syracuse University		H. H. Caspers, Code RRIS	1
Syracuse 10, N. Y.		H. Piller, Code RRIS	1
Attn: Dr. N. Ginsburg, Physics		C. P. Haber, Code RRC	1
Department	1	F. C. Essig, Code RRE	1
H. Levenstein, Solid State		B. F. Husten, Code RF	1
Laboratory	1	H. A. Bulgerin, Code RFV	1
University of Arizona		R. L. Bauer, Code RFS	1
Lunar and Planetary Laboratory		J. J. Nastronero, Code RFA	1
Tucson, Ariz.		F. C. Alpers, Code RM-2	1
Attn: R. B. Goranson	1	T. B. Jackson, Code RMI	1
University of California		A. C. Sevey, Code RMD	1
Riverside, Calif.		Technical Library, Code ADIL	2
Attn: A. Lawson	1		
University of Michigan			
Willow Run Laboratories			
P. O. Box 618			
Ann Arbor, Mich.			
Attn: W. Wolfe	1		
Washington University			
St. Louis, Mo.			
Attn: Dr. E. U. Condon	1		

NOLC 9-63-160

<p>Naval Ordnance Laboratory Corona. (NOLC Report 589) AN INSTRUMENT FOR MEASURING THE SPECTRAL EMITTANCE OF A SOLID IN THE INFRARED FROM 16 TO 45 MICRONS AT LOW TEMPERATURES, by D. L. Stierwalt, D. D. Kirk, and J. B. Bernstein, Research Department. 15 August 1963. 17 pp. UNCLASSIFIED</p> <p>As part of a continuing research project for measuring the emittance of solids, a Beckman IR-3 spectrophotom- eter has been modified. This modification permits the measurement of spectral emittance of solids in the 16 to 45 <math>\mu</math> region in the temperature range between 77 and 473°K.</p>	<p>I. Spectrophotometers - Performance I. Beckman IR-3 spectro- photometer II. Title III. Title: Spectral emit- tance of a solid in the infrared. IV. Stierwalt, D. L. Kirk, D. D. Bernstein, J. B.</p> <p>Task: NASA Purchase Order No. W-11, 400-B</p>	<p>Naval Ordnance Laboratory Corona. (NOLC Report 589) AN INSTRUMENT FOR MEASURING THE SPECTRAL EMITTANCE OF A SOLID IN THE INFRARED FROM 16 TO 45 MICRONS AT LOW TEMPERATURES, by D. L. Stierwalt, D. D. Kirk, and J. B. Bernstein, Research Department. 15 August 1963. 17 pp. UNCLASSIFIED</p> <p>As part of a continuing research project for measuring the emittance of solids, a Beckman IR-3 spectrophotom- eter has been modified. This modification permits the measurement of spectral emittance of solids in the 16 to 45 <math>\mu</math> region in the temperature range between 77 and 473°K.</p>	<p>I. Spectrophotometers - Performance I. Beckman IR-3 spectro- photometer II. Title III. Title: Spectral emit- tance of a solid in the infrared. IV. Stierwalt, D. L. Kirk, D. D. Bernstein, J. B.</p> <p>Task: NASA Purchase Order No. W-11, 400-B</p>	<p>I. Spectrophotometers - Performance I. Beckman IR-3 spectro- photometer II. Title III. Title: Spectral emit- tance of a solid in the infrared. IV. Stierwalt, D. L. Kirk, D. D. Bernstein, J. B.</p> <p>Task: NASA Purchase Order No. W-11, 400-B</p>	<p>I. Spectrophotometers - Performance I. Beckman IR-3 spectro- photometer II. Title III. Title: Spectral emit- tance of a solid in the infrared. IV. Stierwalt, D. L. Kirk, D. D. Bernstein, J. B.</p> <p>Task: NASA Purchase Order No. W-11, 400-B</p>
This card is UNCLASSIFIED		This card is UNCLASSIFIED		This card is UNCLASSIFIED	